

# INTRODUCING....

## The **New IAC 350** Series Maxi Sound Shelter



The 350 Series Maxi Booth is IAC's newest Audiometric Booth. With internal dimensions of 894mm wide x 948mm deep x 2000mm high, the 350 Maxi is an oversized version of IAC's best selling Booth—the 250 Series Mini Sound Shelter.

The 350 Series Maxi Sound Shelter is the ideal environment to carry out hearing screening tests, whether it be at a healthcare practice or annual employee evaluation testing for industry.



### **Standard Features:-**

- Available from stock
- Supplied in 'kit' form
- Fitted with castors or anti-vibration mounts
- Left hand or right hand door option
- BS EN ISO 9001/14001 quality registered
- Jack panel incorporating USB connection
- Maintenance free LED lighting
- Ventilation
- Polyester powder coated RAL 9002 (Grey/White)
- External Dimensions: 1000mm wide x 1069mm deep x 2253mm high

**Contact IAC for more information  
or a no-obligation quotation...**



IAC House, Moorside Road, WINCHESTER,  
Hampshire, S023 7US, United Kingdom  
Tel: +44 (0) 1962 873000 Fax: +44 (0) 1962 873111  
Email: [info@iacl.co.uk](mailto:info@iacl.co.uk)  
Website: [www.industrialacoustics.com](http://www.industrialacoustics.com)